

Antimicrobial Resistance Panel 1: *Acinetobacter baumannii* Lipid A/Fatty Acid Pathway

Catalog No. NR-55640

Product Description:

NR-55640 consists of a 2-member panel of laboratory generated deletion mutant strains of *Acinetobacter baumannii* (*A. baumannii*), strain ATCC® 19606™, which are defective in lipid A biosynthesis or lipopolysaccharide (LPS) transport to the outer cell wall.

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

Table 1: Kit Components

COMPONENT NUMBER	DESCRIPTION	LOT NUMBER	DATE OF MANUFACTURE
NR-51929	<i>A. baumannii</i> , strain NB48062-TMT0028	70043418	07APR2021
NR-51938	<i>A. baumannii</i> , strain NB48062-LMD0007	70043416	14APR2021

NR-51929 and NR-51938 were produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Table 2: *A. baumannii*, Strain NB48062-TMT0028 (NR-51929)

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-negative rods	Gram-negative rods
Colony morphology 2 days at 37°C in an aerobic atmosphere on Tryptic Soy agar	Report results	Circular, convex, entire, smooth and cream (Figure 1a)
Growth at 44°C ± 2°C ¹ 1 day in an aerobic atmosphere on Tryptic Soy agar	Growth	Growth
Motility BBL™ Motility Test Medium w/TTC Indicator for 1 day at 37°C in an aerobic	Report results	Non-motile
VITEK® MS (MALDI-TOF)	<i>A. baumannii</i>	<i>A. baumannii</i> (99.9%)
Antibiotic Susceptibility Profile²		
BD BBL™ Sensi-Disc™ susceptibility test disc 1 day at 35°C in an aerobic atmosphere on Mueller Hinton agar Novobiocin	Report results	30 to 33 mm
Etest® antibiotic test strips 1 day at 35°C in an aerobic atmosphere on Mueller Hinton agar		
Azithromycin	Report results	0.5 to 1.5 µg/mL
Rifampin	Report results	≤ 0.002 µg/mL
Ciprofloxacin	Report results	Sensitive (0.25-0.5 µg/mL)
Fusidic acid	Report results	0.125 µg/mL
Polymyxin B	Report results	Resistant (≥ 1024 µg/mL)
Kanamycin	Report results	Pending

TEST	SPECIFICATIONS	RESULTS
Genotypic Analysis Digital DNA-DNA hybridization (dDDH) ³ Presence of kanamycin resistance cassette (<i>kanR2</i>)	≥ 70% for species identification <i>kanR</i> gene present	<i>A. baumannii</i> (100%) <i>kanR</i> gene present
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Growth at 44°C differentiates *A. baumannii* from *A. calcoaceticus* and *A. pittii*, which do not grow at 44°C.

²Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

³Relatedness between bacterial strains has traditionally been determined using DDH. For additional information, refer to Auch, A. F., et al. "Digital DNA-DNA Hybridization for Microbial Species Delineation by Means of Genome-to-Genome Sequence Comparison." *Stand. Genomic Sci.* 2 (2010): 117-134. PubMed: 21304684.

Table 3: *A. baumannii*, Strain NB48062-LMD0007 (NR-51938)

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Growth at 44°C ± 2°C ¹ 1 day in an aerobic atmosphere on Tryptic Soy agar Motility BBL™ Motility Test Medium w/TTC Indicator for 1 day at 37°C in an aerobic atmosphere VITEK® MS (MALDI-TOF)	Gram-negative rods Report results Growth Report results <i>A. baumannii</i>	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1b) Growth Non-motile <i>A. baumannii</i> (99.9%)
Antibiotic Susceptibility Profile² BD BBL™ Sensi-Disc™ susceptibility test disc 1 day at 35°C in an aerobic atmosphere on Mueller Hinton agar Novobiocin Etest® antibiotic test strips 1 day at 35°C in an aerobic atmosphere on Mueller Hinton agar Azithromycin Rifampin Ciprofloxacin Fusidic acid Polymyxin B Kanamycin	Report results Report results Report results Report results Report results Report results	30 mm 0.75 to 1.5 µg/mL ≤ 0.002 µg/mL Sensitive (0.16 to 0.25 µg/mL) ≤ 0.125 µg/mL Resistant (512 to 768 µg/mL) Pending
Genotypic Analysis Digital DNA-DNA hybridization (dDDH) ³ Presence of kanamycin resistance cassette (<i>kanR2</i>)	≥ 70% for species identification <i>kanR2</i> gene present	<i>A. baumannii</i> (100%) <i>kanR2</i> gene present
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Growth at 44°C differentiates *A. baumannii* from *A. calcoaceticus* and *A. pittii*, which do not grow at 44°C.

²Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

³Relatedness between bacterial strains has traditionally been determined using DDH. For additional information, refer to Auch, A. F., et al. "Digital DNA-DNA Hybridization for Microbial Species Delineation by Means of Genome-to-Genome Sequence Comparison." *Stand. Genomic Sci.* 2 (2010): 117-134. PubMed: 21304684.

Figure 1a: NR-51929
Colony Morphology




Figure 1b: NR-51938
Colony Morphology



/Heather Couch/
Heather Couch

Program Manager or designee, ATCC Federal Solutions

26 JAN 2022

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected by ATCC® and the contributor to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

