

**Acanthamoeba sp., Strain CDC:V540**

**Catalog No. NR-46475**

This reagent is the tangible property of the U.S. Government.

**Product Description:** *Acanthamoeba* sp., strain CDC:V540 is a clinical isolate collected in 2003 from the frozen brain tissue of a 60-year-old female patient in Massachusetts, USA.

**Lot<sup>1,2</sup>: 2186**

**Manufacturing Date: 04NOV2016**

TEST	SPECIFICATIONS	RESULTS
<b>Cellular Morphology</b>	Report results	Adherent and non-adherent, refractile
<b>Genotyping</b> Sequencing of 18S ribosomal RNA gene (~ 430 base pairs)	Consistent with <i>Acanthamoeba</i> sp.	Consistent with <i>Acanthamoeba</i> sp. (genotype T1)
<b>Functional Activity by PCR Amplification<sup>3</sup></b> 18S ribosomal RNA gene (amplicon ASA.S1)	423 to 551 base pair amplicon	~ 450 base pair amplicon
<b>Viable Cell Count by Hemocytometry (pre-freeze)</b>	> 10 <sup>6</sup> cells/mL	6.7 × 10 <sup>6</sup> cells/mL
<b>Viability<sup>4</sup></b>	Growth	Growth
<b>Sterility (21-day incubation)</b> Harpo's HTYE broth <sup>5</sup> , 37°C and 26°C, aerobic Tryptic Soy broth, 37°C and 26°C, aerobic Sabouraud Dextrose broth, 37°C and 26°C, aerobic DMEM with 10% FBS, 37°C, aerobic Sheep Blood agar, 37°C, aerobic Sheep Blood agar, 37°C, anaerobic Thioglycollate broth, 37°C, anaerobic	No growth No growth No growth No growth No growth No growth No growth	No growth No growth No growth No growth No growth No growth No growth

<sup>1</sup>NR-46475 was produced by cultivation of the deposited material in Peptone Yeast Glucose (PYG) medium (ATCC® medium 712) for 8 days at 30°C in an aerobic atmosphere to produce this lot.

<sup>2</sup>Quality control testing completed on post-freeze material unless specified as pre-freeze.

<sup>3</sup>PCR amplification was performed using the JDP1 and JDP2 primer set (JDP1: 5'-GGCCCAGATCGTTTACCGTGAA-3' and JDP2: 5'-TCTCACAAAGCTGCTAGGGAGTCA-3') as described in Schroeder, J. M., et al. "Use of Subgenetic 18S Ribosomal DNA PCR and Sequencing for Genus and Genotype Identification of *Acanthamoeba* from Humans with Keratitis and from Sewage Sludge." *J. Clin. Microbiol.* 39 (2001): 1903-1911. PubMed: 11326011.

<sup>4</sup>Viable cells were observed after 1 day at 30°C in an aerobic atmosphere in PYG medium.

<sup>5</sup>Atlas, Ronald M. *Handbook of Microbiological Media*. 3rd ed. Ed. Lawrence C. Parks. Boca Raton: CRC Press, 2004, p. 798.

**Date:** 28 MAR 2017

**Signature:**



BEI Resources Authentication

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

