

SUPPORTING INFECTIOUS DISEASE RESEARCH

Product Information Sheet for NR-19691

Vibrio cholerae Gateway[®] Clone Set, Recombinant in Escherichia coli, Plate 13

Catalog No. NR-19691

This reagent is the tangible property of the U.S. Government.

For research use only. Not for human use.

Contributor:

Pathogen Functional Genomics Resource Center at the J. Craig Venter Institute

Manufacturer:

BEI Resources

Product Description:

Production in the 96-well format has increased risk of crosscontamination between adjacent wells. Individual clones should be purified (e.g. single colony isolation and purification using good microbiological practices) and sequence-verified prior to use. BEI Resources does not confirm or validate individual mutants provided by the contributor.

The Vibrio cholerae (V. cholerae) Gateway® clone set consists of 46 plates which contain 3813 sequence validated clones from V. cholerae, strain El Tor N16961 cloned in Escherichia coli (E. coli) DH10B-T1 cells. Each open reading frame was constructed in vector pDONR™221 with a native start codon and stop codon. The library was independently cloned and sequence verified by the Harvard Institute of Proteomics. Detailed information about each clone is shown in Table 1.

Information related to the use of Gateway[®] Clones can be obtained from Invitrogen[™]. Recombination was facilitated through an attB substrate (attB-PCR product or a linearized attB expression clone) with an attP substrate (pDONR™221) to create an attL-containing entry clone. The entry clone contains recombinational cloning sites, attL1 and attL2 to facilitate gene transfer into a destination vector, M13 forward and reverse priming sites for sequencing and a kanamycin resistance gene for selection. Please refer to the Invitrogen[™] Gateway[®] Technology Manual for additional details.

Material Provided:

Each inoculated well of the 96-well plate contains approximately 60 μ L of *E. coli* culture (strain DH10B-T1) in Luria Bertani (LB) broth containing 50 μ g/mL kanamycin supplemented with 15% glycerol.

Packaging/Storage:

NR-19691 was packaged aseptically in a 96-well plate. The product is provided frozen and should be stored at -80°C or colder immediately upon arrival. For long-term storage, the vapor phase of a liquid nitrogen freezer is recommended. Freeze-thaw cycles should be avoided.

Growth Conditions:

Media:

LB broth or agar containing 50 µg/mL kanamycin

Incubation:

Temperature: *E. coli*, strain DH10B-T1 clones should be grown at 37°C.

Atmosphere: Aerobic

Propagation:

- Scrape top of frozen well with a pipette tip and streak onto agar plate.
- 2. Incubate the plates at 37°C for 1 day.

Citation:

Acknowledgment for publications should read "The following reagent was obtained through BEI Resources, NIAID, NIH: *Vibrio cholerae* Gateway® Clone Set, Recombinant in *Escherichia coli*, Plate 13, NR-19691."

Biosafety Level: 1

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the following publication: U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health. Biosafety in Microbiological and Biomedical Laboratories. 5th ed. Washington, DC: U.S. Government Printing Office, 2009; see www.cdc.gov/biosafety/publications/bmbl5/index.htm.

Disclaimers:

You are authorized to use this product for research use only. It is not intended for human use.

Use of this product is subject to the terms and conditions of the BEI Resources Material Transfer Agreement (MTA). The MTA is available on our Web site at www.beiresources.org.

While BEI Resources uses reasonable efforts to include accurate and up-to-date information on this product sheet, neither ATCC® nor the U.S. Government make any warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. Neither ATCC® nor the U.S. Government warrants that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, use and disposal. ATCC® and the U.S. Government are not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, the U.S. Government, ATCC®, their suppliers and contributors to BEI Resources are not liable for damages arising from the misidentification or misrepresentation of products.

Use Restrictions:

This material is distributed for internal research, noncommercial purposes only. This material, its product or its derivatives may not be distributed to third parties. Except as performed under a U.S. Government contract, individuals

BEI Resources

www.beiresources.org

E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898



SUPPORTING INFECTIOUS DISEASE RESEARCH

Product Information Sheet for NR-19691

contemplating commercial use of the material, its products or its derivatives must contact the contributor to determine if a

license is required. U.S. Government contractors may need a

References:

 Heidelberg, J. F., et al. "DNA Sequence of both Chromosomes of the Cholera Pathogen Vibrio cholerae." Nature 406 (2000): 477-483. PubMed: 10952301

ATCC[®] is a trademark of the American Type Culture Collection.

Table 1: Vibrio cholerae Gateway® Clones, Plate 13

Clone ID	Well Position	ORF Length	Locus ID	Symbol	Product	Accession Number
174254	A02	459	VC1848	cysS	cysteinyl-tRNA synthetase	NP_231482.1
174442	A03	435	VC1110	serS	seryl-tRNA synthetase	NP_230656.1
174127	B01	N/A	VCA0290	rplT	ribosomal protein L20	N/A
174262	B02	461	VC0560	ffh	signal recognition particle protein	NP_230211.1
175196	B03	107	VC0750		hesB family protein	NP_230399.1
174135	C01	120	VC1801		hypothetical protein	NP_231436.1
174582	C02	209	VC2596	rpIC	ribosomal protein L3	NP_232224.1
174588	D02	236	VC2625	rpe	ribulose-phosphate 3-epimerase	NP_232253.1
174531	E01	144	VC2577	rplO	ribosomal protein L15	NP_232185.1
174685	E02	356	VC2499		conserved hypothetical protein	NP_232128.1
174555	F01	167	VC2579	rpsE	ribosomal protein S5	NP_232207.1
174717	F02	370	VC2036	asd	aspartate-semialdehyde dehydrogenase	NP_231670.1
174564	G01	177	VC2581	rpIF	ribosomal protein L6	NP_232209.1
175022	G02	240	VC0710		conserved hypothetical protein	NP_230359.1
174235	H01	440	VC2242	mesJ	mesJ protein	NP_231873.1
175038	H02	253	VC0708		conserved hypothetical protein	NP_230357.1

BEI Resources
www.beiresources.org

E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898