

Certificate of Analysis for NR-871

Canine Tumor Fibroblast (A-72) Cells, Chemically Inactivated Mock-Infected Cell Control

Catalog No. NR-871

This reagent is the property of the U.S. Government.

Product Description: Cell lysate and supernatant from mock-infected canine tumor fibroblast (A-72) cells. The suspension of cell lysate and supernatant was treated with binary ethyleneimine to simulate virus inactivation.

Lot: 4618586 Manufacturing Date: 08FEB2005

TEST	SPECIFICATIONS	RESULTS
Cell Culture Immunofluorescence Assay ¹ on A-72 Cells	Report results	< 10 fluorescent focus units/mL
Antigen-Capture ELISA ²	Report results	< 10
Sterility (21-day incubation) Harpo's HTYE broth ³ , 37°C and 26°C, aerobic Trypticase soy broth, 37°C and 26°C, aerobic Sabouraud broth, 37°C and 26°C, aerobic Sheep blood agar, 37°C, aerobic Sheep blood agar, 37°C, anaerobic Thioglycollate broth, 37°C, anaerobic DMEM with 10% FBS, 37°C and 5% CO ₂	No growth	No growth
Mycoplasma Contamination Agar and broth culture (14-day incubation at 37°C) DNA Detection by PCR of Test Article nucleic acid	None detected None detected	None detected None detected

¹Welch, S. K. and L. J. Saif. "Monoclonal Antibodies to a Virulent Strain of Transmissible Gastroenteritis Virus: Comparison of Reactivity with Virulent and Attenuated Virus." <u>Arch. Virol.</u> 101 (1988): 221-235. PubMed: 2845894.

Date: 26 JAN 2009 **Signature:** Signautre on File

Title: Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected by ATCC® or the contractor to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

Biodefense and Emerging Infections Research Resources Repository P.O. Box 4137

800-359-7370

²Titer is expressed as the reciprocal of the highest dilution that resulted in a mean absorbance greater than the mean absorbance of the mock-infected control plus three standard deviations.

³Atlas, Ronald M. Handbook of Microbiological Media. 3rd ed. Ed. Lawrence C. Parks. Boca Raton: CRC Press, 2004, p. 798.